


Veterinary Scales

Scales for Veterinarians, Pet Stores, Kennels or Laboratories.

The Digital Weight Indicator "Locks-On" Weight for Easy Handling. Battery or AC Power.



Rubber mat shown optional.

VET50 Digital Veterinary Scale

Ideal for weighing cats and medium to small size dogs

50 lb x 0.02 lb / 20 kg x 0.01 kg



VET400 Digital Veterinary Scale

Ideal for weighing large dogs or heavy animals

400 lb x 0.2 lb / 200 kg x 0.01 kg

Detecto Veterinary Scales

Detecto's Vet Scales Utilize the State-of-the-Art Model 750 Digital Weight Indicator

- Large 0.75" Digits - Easy-To-Read LCD.
- LOCK/RELEASE Button Retains Animal's Weight For Efficient Handling.
- Easily-Removable Battery Cover.
- 6 "AA" Batteries or Optional AC Power.
- Durable Metal Construction.
- Mounting Brackets For Desk or Wall.



VET50 Digital Veterinary Scale

This portable, compact-sized scale with a remote weight display is the ideal scale for small animals. Its scale base can be positioned up to 8 feet away from the display. The scale's medium base size allows it to be used on the floor or tabletop and its light weight makes it ultra portable. An optional rubber mat (as pictured below) is also available.

Capacity: 50 lb x 0.02 lb / 20 kg x 0.01 kg
Platform Size: 3" H x 14" W x 12" D / 8 cm H x 36 cm W x 30 cm D
Shipping Weight: 15 lb / 7 kg



VET400 Digital Veterinary Scale

The perfect scale for veterinarians, kennels, labs or anyone handling medium to large size animals. The large, low-profile platform makes it easy for animals to step onto the floor scale without hesitancy, and an optional rubber mat is also available. The removable stainless steel platform makes clean-up fast and easy. The indicator comes with a 8' coiled cable and adjustable mounting bracket for either desk or wall mount convenience.

Capacity: 400 lb x 0.2 lb / 200 kg x 0.1 kg
Platform Size: 1.75" H x 42" W x 22" D / 44 cm H x 107 cm W x 56 cm D
Shipping Weight: 76 lb / 34 kg



The scale and display are battery powered with an automatic sleep mode to conserve battery power when not in use.



Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By:



Medical Measurements
6404 Johnston Street
Lafayette, LA 70503
800-706-9270