

M24

Medical Measurements
Lafayette, LA
800-706-9270
www.medicalmeasurements.com

The M24 centrifuge is a truly affordable microhematocrit centrifuge. This unit is designed for its ease of use in the laboratory. Spin hematocrits in less than 5 minutes at 12,000 rpm. The M24 has a microhematocrit reader disk provided with every unit. Read the tubes directly on the rotor without picking them up.

ADVANTAGES

Fast results in just five minutes
Spins 24 standard 75mm capillary tubes
Ready-to-use...no add-on costs
Quiet, vibration-free operation
Reader disk included

FEATURES

Impact-resistant safe-view lid
All-metal body construction
Cooling fan for motor compartment
Aluminum rotor for low heat transfer
Heavy duty high-speed motor
Safety shut-off switch in lid
Locking lid
Wide mouth for easy clean up

SPECIFICATIONS

12,000 rpm fixed speed
Nominal RCF=14,800 g
Brushed motor
110V AC (220V AC available)

CONTROLS

Illuminated power switch
Electric brake for quick stop
Auto-off 15-minute timer with bell

INCLUDES

24-place hematocrit rotor
Reader disk— read tubes directly on the rotor



M24

Model #	Description
HEC-24HF-7501	Microhematocrit centrifuge

- 24-Place Microhematocrit Centrifuge
- Affordable
- Results in Less Than 5 Minutes

WEIGHT AND DIMS

Height: 10" (235mm)
Length: 13" (330mm)
Width: 12" (304mm)
Weight: 23 lbs. (10.4kg)



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com