

Combo V24

To Purchase or for more information:

Medical Measurements

800-706-9270

www.medicalmeasurements.com

Accomplish both test tube and microhematocrit separations with the new Combo centrifuge. Simply choose a box (with or without a digital tachometer) and choose your rotors. The Combo spins either test tubes, capillary tubes or microtubes by simply changing rotors — no tools required! Achieve hematocrit separations in less than five minutes at 12,000 rpm; then accurately determine results with the included reader disc. Other features include variable speed control, 30-minute timer, rubber feet, and all-metal body construction.

ADVANTAGES

- One unit for all centrifuge needs
- Fast spin, fast stop, fast read
- Variable speed
- Optional tachometer
- Quiet vibration-free operation

FEATURES

- Suction-cupped feet to prevent slipping
- Impact-resistant safe-view lid
- All-metal body construction
- Cooling fan for motor compartment
- Heavy duty high-speed motor
- Safety shut off switch in lid
- Locking lid
- Wide mouth for easy clean-up

SPECIFICATIONS

12,000rpm on High speed setting

- Microhematocrit rotor: 14,500g
- Microtube rotor: 12,900g

4,000rpm on Low speed setting

- Test tube rotor: 2,415g

CONTROLS

- Illuminated power switch
- Electric brake for quick stop
- High/low speed setting switch
- Auto-off 30 minute timer with bell
- Variable speed control

WEIGHT AND DIMS

Height: 11" (277mm)
Length: 17" (432mm)
Width: 14" (355mm)
Weight: 31.3 lbs (14.2kg)



Combo V24

Model #	Description
CMP-BX7N-7701	Combo centrifuge box with variable speed dial, timer, and brake (no tachometer)
CMP-BX7N-77T1	Combo centrifuge box with digital tachometer, variable speed dial, timer, and brake

- Spin Hematocrit or Test Tubes
- Easy-Change Rotors
- Tachometer Option for Exact Forces

CHOOSE THE ROTORS



CMP-RT24-75HE

Spin standard 75mm capillary tubes with the 24-place microhematocrit rotor. Results can be read directly from the rotor with the included reader disk.



CMP-RT24-02MT

The 24-place microtube rotor holds 1.5 to 2 ml sized tubes.



CMP-RT12-15TT

The 12-place test-tube rotor holds 3ml-15ml sized tubes with included tube shields and inserts for smaller tubes. 24 place rotor also available.



CMP-RT24-15TT

The 24-place test-tube rotor holds 15ml sized tubes with included tube shields and inserts for smaller tubes. 12 place rotor also available.



P 770.270.1394

F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com