

T-Edge 10" Autoclave

The only gravity autoclave
upgradable to pre/post vacuum
by software!

- Larger chamber capacity and shorter cycle times
- Convenient and easy loading with door opening to 110°
- Advanced documentation and traceability options and barcode

32 MIN.
WRAPPED
CYCLE TIME
INCLUDING
DRYING!



Remote process
monitoring on
smart devices.


Designed to meet all current
sterilization standards ANSI/AAMI ST55


Tuttnauer
Innovation · Legacy · Partnership

Advanced barcode traceability



Receive barcodes with sterilization data (including expiration date).



Apply barcode label to each pouch.



- Optimal chamber capacity (23 Liter/6 Gallons)
- Convenient and easy loading with door opening to 110°
- Assorted water filling and drainage methods (top, front or automatic)
- Range of automatic one-button programs
- User-friendly touchpad and control panel
- Advanced documentation and traceability options and barcode
- Remote Wi-Fi monitoring via smartphone or tablet

Technical Specifications

Autoclave Model	Chamber Size (In)	Volume	Overall Dimensions D x W x H (In)	Number of trays	Connections	Standard Cassette Capacity	Item Supplies	Power Supply	Power Consumption	Weight
T-Edge 10	18" depth x 10" diameter	6 gal/ 23 liters	22.8"x18.9" x19.23"	5	1 ethernet port, 4 USB ports, WiFi, integrated automatic water filling and drain connection	3 Full + 3 Half	Drain tube, tray holder	230v, Frq 60Hz	2000W	110 lbs.

Sterilization Programs

CLASS S

Icon	Cycle Program	Load Type	Cycle Time with Drying (minutes)
	Unwrapped Instruments 270° F	Unwrapped Instruments	20 Minutes
	Wrapped Pouches 270° F	Wrapped Pouches	32 Minutes
	Unwrapped delicate 250° F	Unwrapped Delicate Instruments	45 Minutes
	Handpieces 270° F	Handpieces	45 Minutes


CLASS B

Icon	Cycle Program	Load Type	Cycle Time with Drying (minutes)
	Unwrapped Instruments 273° F	Unwrapped Instruments	20 Minutes
	Wrapped Pouches 273° F	Wrapped Pouches	39 Minutes
	Unwrapped Delicate 250° F	Unwrapped Delicate Instruments	35 Minutes
	Wrapped Delicate 250° F	Wrapped Delicate Instruments	60 Minutes



Tuttnauer
Innovation · Legacy · Partnership

Learn from our experts.
Visit our blog: tuttnauer.com/blog



Tuttnauer USA Co. Ltd.
25 Power Drive, Hauppauge, NY 11788
Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720
E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

